



innovative
Sensor
Technology

SFS01 Evaluation Kit



Evaluation Kit for gas flow monitoring with a fast calorimetric flow sensor

Pre-calibrated gas flow module with a 1x1 mm PEEK flow channel and choice of digital (I2C) and analog signal output

Product Name: SFS01 EVAKIT

Gehäuseabmessungen: 88.0 x 27.5 x 30.4 mm

Gehäusematerial: PEEK

Anschluss: QSM M5-4

Strömungskanal innerer $\varnothing$: 1x1 mm

Strömungskanal Länge: 35 mm

Wetted materials: PEEK, silicon based chip

Ansprechzeit T63: < 10 ms

Gewicht: 37 g

Speisespannung: 5V DC

Ausgangssignal: Mass flow in sccm, raw signal

Genauigkeit: < 3 % f.s.

Reproduzierbarkeit: < 0.3 % f.s.

Product code: 105059

Product details

SFS01 Evaluation Kit > Details

The online shop

Quantity (pieces)	Price (per piece)
--------------------------	--------------------------

1-5	CHF 339.90
-----	------------

Lager: **13**