



innovative
Sensor
Technology

SFS01 Evaluation Kit



Evaluation Kit for gas flow monitoring with a fast calorimetric flow sensor

Pre-calibrated gas flow module with a 1x1 mm PEEK flow channel and choice of digital (I2C) and analog signal output

Product Name: SFS01 EVAKIT

Housing dimensions: 88.0 x 27.5 x 30.4 mm

Housing material: PEEK

Fluidic connection: QSM M5-4

Flow channel inner dimensions: 1x1 mm

Flow channel length: 35 mm

Wetted materials: PEEK, silicon based chip

Response time t_{63} : < 10 ms

Weight: 37 g

Supply voltage: 5V DC

Output signal: Mass flow in sccm, raw signal

Accuracy: < 3 % f.s.

Repeatability: < 0.3 % f.s.

Product code: 105059

Product details

SFS01 Evaluation Kit > Details

The online shop

Quantity (pieces)	Price (per piece)
--------------------------	--------------------------

1-5	CHF 339.90
-----	------------

Stock: 13