



innovative  
Sensor  
Technology

## Thermal gas flow sensor FS7



Flow range from 0 to 100 m/s, -20 °C to +150 °C

For measuring flow rates using a symmetrical heater design and having an excellent sensitivity and reproducibility. This all-rounder features long-term stability and reproducibility.

**Product Name:** FS7.0.1L.195

**Measured medium:** Gas Flow

**Flow direction:** One direction (anemometric)

**Operating temperature range:** -20 °C to +150 °C

**Chip size/dimensions:** 6.9 x 2.4 x 0.20 mm

**Connection type:** Cu/Ag stranded wire, PTFE-insulated

**Wire dimensions:** 195 mm

**Flow range CTA mode:** 0/ms to 100 m/s

**Response time t63:** ~200 ms (jump from 0 to 10000 sccm)

**Response sensitivity:** 0.01 m/s

**Substrate:** Ceramic

**Sensor wire:** 3 pins, Cu/Ag-stranded wires, PTFE insulated

**Product Old code:** 050.00216

**Product code:** 103705

## Product details

### Thermal gas flow sensor FS7 > Details

#### Product description FS7

The IST AG FS7 thermal mass flow sensors are the successors of the FS5 flow sensors and exhibit a symmetrical heater design and enhanced sensitivity. The FS7 flow sensor is applicable in gas and offer excellent long-term stability. The small thermal mass of the FS7 sensor provides a fast response time. Additionally, the standard FS7 sensor is also available with a housing which is easily implemented into various applications.

The IST AG FS7 flow sensors are optimal for measuring gas flow in various flow applications.

#### Benefits & Characteristics of the FS7 gas flow sensor

- Easy adaptation in various applications and housings
- Excellent long-term stability
- Simple signal processing
- Simple calibration
- Stable platinum technology
- Excellent reproducibility

[Customer-specific sensor available upon request.](#)

#### The online shop

<b>Quantity (pieces)</b>	<b>Price (per piece)</b>
--------------------------	--------------------------

5-49	CHF 28.60
------	-----------

50-99	CHF 22.70
-------	-----------

100-150	CHF 20.09
---------	-----------

Stock: **134**