



innovative
Sensor
Technology

Biosensor LV5 , packaged sensor



Flow-through biosensor for multiparametric measurement (Glucose and Lactate as standard).

Note: Not for medical, diagnostics and use on humans. For evaluation use only. For more information contact IST AG.

Product Name: B.LV5.GL.C010.ST.L050.ECC.R

Analyte: Glucose, Lactate

Sensor dimensions: 10.0 x 4.0 x 0.4 mm (H with flow cell: 1.23 mm)

Construction type: Flow-through cell

Flow range rating: 0.3 μ l/min to 10 ml/min

Fluidic connection: \varnothing inner: 0.5 mm (0.15 mm available upon request)

Working electrode: Platinum covered with enzyme membrane

Blank electrode: For background compensation

Reference electrode: Silver / silver chloride

Counter electrode: Platinum

Glucose measurement range: 0.05 mM to 25 mM / 0.01 to 4.5 g/l

Lactate measurement range: 0.02 mM to 15 mM / 0.002 to 1.5 g/l

Sterilization: Beta or gamma sterilization is possible (25 kGy typical); sensitivity increases upon radiation; lactate sensors life time is reduced; Initial bioburden < 1 cfu per sensor

Shelf life: >6 months (at recommended storage conditions)

Operational life time: > 2 weeks continuous operation; > 4 weeks in analyzer mode

Product code: 105128

Product details

Biosensor IV4 for Glucose measurement > Details

Benefits & Characteristics of the LV5 sensor

- Excellent long-term stability over a month in continuous measurement mode (lifetime varies depending on various factors, including the analyte concentration and buffer system (medium) amongst others)
- Long shelf life (under appropriate environmental conditions)
- Gamma and beta sterilization compatible
- Suitable for continuous and analyzer measurement modes
- Fast response time (depends on the buffer system (medium))
- Reference, counter and blank electrodes on-chip
- Integrated flow cell (various volumes)

Disclaimer

Evaluation product for professionals to be used solely for research and development purposes! Not for medical and diagnostic use. Not to be used on humans. For more information contact IST AG.

Quality

Consistent with the well-known, high-quality standards in Switzerland, IST AG is certified according to ISO 9001:2015 (quality) and ISO 14001:2015 (environment).

Appropriate processes are part of our daily work. They are regularly audited and extended parallel to the growth of our company.

[> Read more](#)

The online shop

Quantity (pieces)	Price (per piece)
--------------------------	--------------------------

1-9	CHF 357.70
-----	------------

10-49	CHF 316.50
-------	------------

50-100	CHF 277.80
--------	------------

Stock: **0**