

physical. chemical. biological.



TSic 501F



Temperature Sensor IC









Benefits & Characteristics

- Fully calibrated
- Outstanding accuracy of ±0.1 K
- Very low power consumption
- Excellent long-term stability
- Accuracy range of 40 K can be shifted (default: +5 °C to +45 °C)

Illustration

L

L2

Technical Data

Dimensions (L / L2 x W x H in mm): 2)	17.30 / 3.81 x 4.57 x 2.3
Operating temperature range:*	-10 °C to +60 °C (-7 °C to +57 °C guaranteed)
Accuracy:*	±0.1 K in the range of +5 °C to +45 °C
Resolution:*	0.034 K
Sampling rate:*	10 Hz
Supply voltage:	$V_{dd} = 3 \ V$ to 5.5 V, high precision operation in range $V_{dd} = 4.5 \ V$ to 5.5 V
Supply current:	typ. 30 μA at 25 °C and $V_{dd} = 3.3 \text{ V}$ for minimal self-heating
Packaging:*	TO92
Output signal:	Analog - see application note ATTSic_E

²⁾ For tolerances, see Application Note



physical. chemical. biological.












Product Photo



Pin Assignment



	Pin 1	Pin 2	Pin 3
TO92	GND	Signal	$V_{dd'}$ Supply voltage (3 V to 5.5 V)

Absolute maximal ratings

	Min	Max
Supply voltage (V _{dd})	-0.3 V	6 V
Voltages to analog I/O – Pins (V_{SIG}, V_{GND})	-0.3 V	V_{dd} +0.3 V
Storage temperature range (T _{STOR})	-10 °C	+60 °C
Non-operating temperature range		

Operating conditions

	Min	Тур	Max
Supply voltage to GND (V+)	2.97 V	5 V	5.5 V
Supply current (I_{Vdd}) at V_{dd} = 3.3 V, RT	25 μΑ	30 μΑ	60 μΑ
Operating temperature range (T _{amb})	-10 °C		+60 °C
Output load capacitance (C _L)			15 nF
External capacitance between V_{dd} and $GND^{1)}$	100 nF (recomme	nded)	
Output load resistance between signal and GND (or $V_{\rm dd}$)	47 kΩ		

 $^{^{\}mbox{\tiny 1)}}\mbox{Recommended}$ as close to TSic $\mbox{V}_{\mbox{\tiny dd}}$ and GND-Pins as possible



physical. chemical. biological.











Temperature accuracies²⁾

T1: +5 °C to +45 °C	±0.1 K
T2: -10 °C to +60 °C	±0.2 K

 $^{2)}$ The sensor is calibrated at 5 V. The provided accuracy is applicable for a supply voltage between 4.5 V and 5.5 V. The accuracy is smaller with a supply voltage between 2.97 V and 4.5 V. For applications where the best accuracy at 3 V is requested, ask for a custom specific, 3 V calibrated device. Other TSic products with custom specific calibrations are available upon request e.g. other temperature range for high accuracy. Accuracy at delivery; the assembly method can influence the accuracy!

Order Information

Description:	Item number:	Former main reference:
TSic 501 TO92	103491	030.00046

Additional Documents

	Document name:
Application Note:	ATTSIC_E

