

physical. chemical. biological.











P0K1.281.2K.B.150.R.S


Platinum sensor with round ceramic housing For low temperatures

Benefits & Characteristics

- Same dimensions as a traditional wire wound sensor easy interchangeability into existing applications
- Easy to assemble (boreholes etc.)
- Excellent long-term stability

- Vibration and temperature shock resistant (depending on assembly)
- Fast response time (depending on assembly)
- Low self-heating

Illustration



Technical Data

Operating temperature range:	-50 °C to +200 °C
Nominal resistance:*	100 Ω at 0 °C
Characteristics curve:*	3850 ppm/K
Long-term stability:	< 0.04 % at 1000 h at maximal operating temperature
Tolerance class (dependent on temperature range):*	IEC 60751 F0.3 B (IST AG reference)
Connection:*	Cu/Ag Stranded wire, AWG26, PTFE insulated, 150 mm long, ends 4 mm stripped
Recommended applied current:2)	1 mA at 100 Ω
²⁾ Self-heating must be considered	
Dimensions:	2.8 x 13.0
Tolerance (chip):	Ø ± 0.2 mm, L ± 1 mm, L $_{\rm W}$ (up to 30 mm) ± 1 mm

Product Photo





Former main reference:

010.00796

physical. chemical. biological.













Order Information

Description:

Application Note:

Document name:

Item number:

100481

ATP_E



