



innovative
Sensor
Technology

SFS01 Evaluation Kit



Evaluation Kit for gas flow monitoring with a fast calorimetric flow sensor

Pre-calibrated gas flow module with a 1x1 mm PEEK flow channel and choice of digital (I2C) and analog signal output

Product Name: SFS01 EVAKIT

□□□□ : 88.0 x 27.5 x 30.4 mm

□□□□ : PEEK

□□□□ : QSM M5-4

□□□□□ : 1x1 mm

□□□□ : 35 mm

Wetted materials: PEEK, silicon based chip

□□□□ **T63:** < 10 ms

□□ : 37 g

□□□□ : 5V DC

□□□□ : Mass flow in sccm, raw signal

□□ : < 3 % f.s.

□□□ : < 0.3 % f.s.

Product code: 105059

Product details

SFS01 Evaluation Kit > Details

The online shop

Quantity (pieces) Price (per piece)

1-5 € 360.41

 : **13**