



innovative
Sensor
Technology

SFS Module 1000sccm EvalKit



Thermal mass flow module with linearized digital / analog output

Fast and compact flow module, suitable for gas dosing systems, manifolds and analysers; Practical Evaluation Kit with fluidic block adapter for direct tube attachment, M5-4 connectors for 4mm tubing included (exchangeable)

Product Name: SFS01-M-01000

□□□□□ : +/-1000sccm (ml/min)

□□□ : Gas Flow

□□□□□ : 0 to 80 °C

□□□□ : Aluminum, FR4

□□□□ : M5-4mm tube adapter or Ø 2.2 mm i.D. flat gasket (w/o fluidic block)

□□□□□ : 2.2 mm

Wetted materials: Aluminum, gold, FR4, silicon, silicon nitride, epoxy, nitrile, brass

□□□□ : PCB connector JST SM08B-SRSS-TB

□□□□ **T63:** 20 ms

□□□□ : I2C or 0 to 1 V

□□ : 2% f.s. (digital output)

 : 0.2 % f.s.

Product code: 156602

Product details

SFS Module 1000sccm EvalKit > Details

The online shop

Quantity (pieces)	Price (per piece)
-------------------	-------------------

1-5	CHF 278.10
-----	------------

 : 0